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Professional Instruction:
ARTZT vitality® Kinesiological Tape back

1. Application

Thanks to their adhesive which is applied in a wavy pattern, the tapes adhere pleasantly to the skin. The tapes can and should be worn for several days. After approximately 5-7 days, their position should be examined and the tapes should be reapplied if necessary. More extended wear has no harmful consequences, even without examination.

The tapes can be replaced as often as desired. Verifying that the skin is sound and unbroken is recommended prior to each application.

2. Areas of application

Back (spine).

In case of radiating pain, especially in combination with sensitivity or motor disturbances, consulting a doctor or physiotherapist is highly recommended.

3. Precautions, safety information

When the complaint is first noticed, the cause should be examined medically in order to exclude the possibility of causes that require medical treatment.

The attending doctor or a physiotherapist should be consulted in case of severe or long-term complaints.

4. Particularities of the tapes

The tapes do not contain harmful additives. The adhesive that is used has been dermatologist tested and is considered especially skin-friendly. This is why there is also normally no risk for persons allergic to conventional plasters. Combinations with other types of therapy are possible and/or recommended (supporting therapy).

Thanks to the high level of wearer comfort of the tapes, all day-to-day work and sports activities can and/or should be conducted as usual with the tapes in place. The same applies to showering. The material is breathable, permeable to moisture and liquids and dries quickly.

5. Application method and removal of the tapes

The skin must be dry and free of grease. Hair should be trimmed (long-hair trimmer) or gently removed entirely by shaving (if possible, the day before applying the tapes). Intensively pressing the tapes in place and rubbing (especially along the edges) after application ensures the adhesive will adhere properly.

To remove the tapes, carefully and gradually pull in the direction of hair growth. The tapes are especially easy to remove in the shower or while taking a bath. Skin cream or body oil help relax the skin after removing the tapes.

The back tapes consist of two strips approximately 40 cm in length. However, the length is individually adjustable. After adjusting the length, the edges have to be rounded again. Both strips are applied to the back depending on the symptoms – that is, according to the area of the complaint.

Because of this positioning, you need a helper in order to apply the tapes.

Remove the two strips from the package and assume an active, upright position which allows you to complete the following activities unencumbered and without pain: Head and arm movements, swallowing and breathing.

1) Remove the backing paper from the first tape strip. Avoid touching the adhesive surface during this process. The two strips are applied to the back so they cross over each other. In order to do so, proceed as follows:

2) Apply the tape to the right shoulder (centred) so that it is slightly offset towards the front and secure this starting point.

3) Now stretch the tape and apply it diagonally (towards the left hip) with tension across the back to the spine by always pressing down the section you have just applied. Only the end of the tape is applied without tension, but also has to be pressed down firmly. Then press the tape down over its entire length. Make sure no folds or ripples form in the process.

4) Carefully remove the backing paper from the other tape strip. Avoid touching the adhesive surface during this process. Apply the tape strip to the left shoulder (centred) so that it is slightly offset towards the front and secure this starting point.

5) Apply the tape on the left shoulder (centered) so that it is slightly offset towards the front and secure this starting point as described under point 2.

6) Similar to the first tape strip, stretch the tape and apply it diagonally (towards the right hip) across the back to the spine as described under point 2.

6. Precautions, safety information

- » Always cut the tapes to size prior to application, never while you are applying them.
- » Always be careful when pulling clothing over the edges of the tapes; avoid loosening the tapes.
- » Avoid the formation of folds and ripples, especially where they could cause pressure.
- » If your complaints persist or worsen, consult your doctor and / or physiotherapist. Never apply the tapes if your complaints have not been diagnosed.
- » Remove the tapes immediately in case of intolerance.
- » Store the material at room temperature and avoid direct exposure to sunlight.

7. Side effects

There are no known side effects with correct application.

Note:

Treatment with kinesiology tape is based on methods from training theory and kinesiology studies, the effects of which are not proven by means of conventional medical studies but rather are assumed based on experience reports that are not evidence-based as well as application observations.

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ARTZT vitality® is a registered trademark of Ludwig Artzt GmbH. Unauthorized use is strictly prohibited.
© 2017 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.

Authors:
Jürgen Förster,
Lehrphysiotherapeut
Christian Röhrs,
Sportphysiotherapeut DOSB
Astrid Buscher,
Dipl. Sportökonomin